

S17 Table

Effects of Retesting Urine after Boiling for ImmuView-BinaxNOW Discordant Result

Specimens for Upenn and SSI Sites^a

Result Change After Boiling ^b	<i>S. pneumoniae</i>		<i>L. pneumophila</i>	
	ImmuView	BinaxNOW	ImmuView	BinaxNOW
FN to FN	2	0	1	4
FN to TP	1	1	1	1
TP to FN	0	0	3	0
TP to TP	1	1	7	2
TN to TN	3	3	0	0
TN to FP	3	0	0	0
FP to FP	2	1	0	0
FP to TN	4	0	1	0

^aThese are summary results and not paired comparisons. The UPenn site tested five *L. pneumophila* and six *S. pneumoniae* ImmuView-BinaxNOW discordant specimens with both BinaxNOW and ImmuView tests. The SSI site tested eight *L. pneumophila* ImmuView-BinaxNOW discordant specimens using only the ImmuView test for seven urines because of limited urine volumes and both BinaxNOW and Immuview for two specimens, and eight *S. pneumoniae* discordant specimens using only the Immuview test because of limited urine volumes; ^b FN, false-negative; TP, true-positive; TN, true-negative; FP, false-positive; FN to FN means false-negative before and after boiling, FN to TP means false-negative before boiling and true-positive after boiling, etc.